Morgan® Lens Instructional Chart

[DA] - Instruktionsdiagram til Morgan® Lens [AR] المخطط الإرشادي لاستخدام عدسة Morgan® Lens [DE] - Morgan® Lens Anwendungstabelle [EL] - Διάγραμμα οδηγιών Morgan® Lens [ES] - Tabla de instrucciones sobre el uso de Morgan® Lens [ET] - Morgan® Lens juhend [FI] - Morgan® Lens Ohjekortti [FR] - Tableau du mode d'emploi Morgan® Lens [HU] - Morgan® Lens útmutató diagram [IT] - Morgan® Lens istruzioni [LT] - „Morgan® Lens“ instrukcijų lentelė [NL] - Morgan® Lens Instructiediagram [NO] - Morgan® Lens-instruksjonskart [PL] - Karta instruktażowa Morgan® Lens [SK] - Inštruktážna tabuľka pre šošovky Morgan® Lens [SV] - Anvisningsblad för Morgan® Lens [TR] - Morgan® Lens Talimat Tablosu [ZH] - Morgan® Lens 指示圖

Contents

[EN] - Morgan® Lens Instructional Chart – for emergency hands-free ocular irrigation	2
[AR] لغسل العين في حالة الطوارئ دون استعمال اليدين Morgan® Lens المخطط الإرشادي لاستخدام عدسة	5
[DA] - Instruktionsdiagram til Morgan® Lens – til nødsituationer med håndfri øjenskyllning.....	8
[DE] - Morgan® Lens Anwendungstabelle – Freihändig nutzbares Augenspülsystem für Notfälle.....	11
[EL] - Διάγραμμα οδηγιών Morgan® Lens - για οφθαλμική έκπλυση έκτακτης ανάγκης χωρίς χέρια.....	14
[ES] - Tabla de instrucciones sobre el uso de Morgan® Lens para la irrigación ocular de emergencia «manos libres»	17
[ET] - Morgan® Lens juhend silmade käed-vabad niisutamiseks hädaolukordades.....	20
[FI] - Morgan® Lens Ohjekortti – hätätilanteessa tapahtuvaan silmän huuhteluun kädet vapaana	23
[FR] - Tableau du mode d'emploi du dispositif d'irrigation oculaire d'urgence mains-libres Morgan® Lens.....	26
[HU] - Morgan® Lens útmutató diagram – Sürgősségi érintés nélküli szemöblítés	29
[IT] - Morgan® Lens istruzioni per l'irrigazione oculare d'emergenza a mani libere	32
[LT] - „Morgan® Lens“ instrukcijų lentelė – skirta skubiam akių drėkinimui nenaudojant rankų	35
[NL] - Morgan® Lens Instructiediagram - voor handenvrije oculaire irrigatie in noodgevallen	38
[NO] - Morgan® Lens-instruksjonskart – for nødsituasjoner med håndfri skylling av øyet	41
[PL] - Karta instruktażowa Morgan® Lens – do doraźnej irygacji oczu bez użycia rąk.....	44
[SK] - Inštruktážna tabuľka pre šošovky Morgan® Lens – na núdzový výplach očí bez použitia rúk.....	47
[SV] - Anvisningsblad för Morgan® Lens – för handsfree ögonspolning i nödsituationer.....	50
[TR] - Morgan® Lens Talimat Tablosu – Acil durumlarda eller kullanılmadan oküler irrigasyon için.....	53
[ZH] - Morgan® Lens 指示圖 – 用於緊急免持式眼灌洗	56

MorTan, Inc.
329 East Pine Street.
Missoula, MT 59802 USA
Tel: 1-800-423-8659 Fax: 1-406-728-9332
mortan@morganlens.com

0413

morganlens.com

[EN] - Morgan® Lens Instructional Chart – for emergency hands-free ocular irrigation

Instructions for using the Morgan Lens for continuous irrigation or delivery of medication to the cornea and conjunctiva.

The Morgan Lens, Morgan Lens Delivery Set, Medi-Duct and MorTan's packaging materials are not made with natural rubber latex.

1 

INSERTION
Instill topical ocular anesthetic, **if available**.

2 

Peel open the sterile packaging and attach the Morgan Lens® to the Morgan Lens Delivery Set® (or a fluid administration set or syringe).

3 

Using solution of choice*, **START A MINIMAL FLOW**. This allows Lens to “float” over cornea and sclera.

4 

Have patient look down, insert Morgan Lens under upper lid.
Have patient look up, retract lower lid, drop Lens in place

5 

Adjust flow to the desired rate and absorb outflow with the Medi-Duct®. Tape tubing to the patient's forehead to prevent accidental removal of the Lens.
DO NOT RUN DRY.

6 

REMOVAL – CONTINUE FLOW.
Have patient look up, retract lower lid — hold position and slide the Lens out. Once Lens is removed, **stop the flow of solution.**

Experts recommend continuous irrigation until the pH returns to normal – approximately 7.0-7.3. After 10 minutes, recheck pH and repeat process if necessary. When irrigation is complete, dispose of Lens properly.

Contraindications: Do not use if there is a protruding foreign body, a penetrating eye injury, or with a suspected or actual rupture of the globe.

Precautions:

- Rx Only – Federal (USA) law restricts this device to sale by or on the order of a physician.
- Device should only be used by professionals trained in ocular irrigation procedures.
- Contents supplied sterile using irradiation. Do not use if package is opened or damaged.
- Single use only. Do not re-use or re-sterilize. Reuse of the device could result in bacterial or viral infection. It could also result in deterioration or damage to the device.
- Dispose after use as biohazard per applicable law and standard medical practice

Notice: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the regulatory authority of the region in which the user and/or patient is established.

Morgan Lens Uses	Solution	Mode with Morgan Lens	Rate	Frequency
Acid burns or other ocular injury due to solvents, gasoline, detergents, etc.	Lactated Ringer's* Solution Recommended due to pH: LR 6.0 to 7.5	Morgan Lens Delivery Set or I.V. set-up	500 ml rapid/free flow. Reassess and continue at slower rate.	Once. Repeat as necessary. Continuous until pH of cul-de-sac is returned to normal (approximately 7.0 to 7.3). Repeat as necessary.
Alkali burns	NS 4.5 to 7.0	Morgan Lens Delivery Set or I.V. set-up	2000 ml rapid/free flow. Reassess. Continue at 50 ml/hour or 15 drops/minute.	Once. Repeat as necessary.
Non-embedded foreign bodies		Morgan Lens Delivery Set or I.V. set-up	500 ml rapid/free flow. Reassess and continue at slower rate.	Once. Repeat as necessary.
Foreign body sensation with no visible foreign body	20 cc sterile solution	20 cc syringe	Slowly without force.	Once. Repeat once if necessary.
Routine pre-operative	10 cc of preferred ocular antiseptic	10 cc syringe	Slowly without force.	Once.
Eyelid surgery (protecting the cornea during eyelid surgery)	Lactated Ringer's* Solution	Morgan Lens Delivery Set or I.V. set-up	4 drops/minute.	During entire procedure.
Severe infection	Lactated Ringer's* Solution with suitable antibiotic and steroid**	Morgan Lens Delivery Set or I.V. set-up	50 ml/hour or 15 drops/minute.	Continuous for 70 hours, then 10- hour intervals until marked improvement.

*MorTan recommends lactated Ringer's because its pH level of 6.0 to 7.5 is much closer to tears (approximately 7.1) than Normal Saline (4.5 to 7.0). Also, the lactate ion in lactated Ringer's exhibits a buffering capacity, returning the pH of either an acid or a base to neutral much more rapidly than a solution (such as saline) without any buffering capacity. **Use with antibiotics and steroids only when indicated.

Video and PowerPoint available online, or call to order a free copy. To place an order or to inquire about our training materials, call us toll free at 1-800-423-8659.

CALL FOR QUANTITY DISCOUNTS

morganlens.com

0413

The Morgan Lens
ZiMorTan Inc.

MorTan, Inc.
329 East Pine Street.
Missoula, MT 59802 USA
Tel: 1-800-423-8659 Fax: 1-406-728-9332
TO ORDER: 1-800-423-8659
email: mortan@morganlens.com

The Morgan Lens, Morgan Lens Delivery Set and the Medi-Duct combine to provide the most effective, hands-free ocular irrigation system available.

The Morgan Lens®


provides hands-free ocular irrigation and/or medication to the cornea and conjunctiva. The Morgan Lens may be used for hours, or even days, providing the continuous treatment that is often required for serious chemical burns.

The Morgan Lens Delivery Set®

allows the simultaneous irrigation of both eyes. Instead of wasting valuable time and money on separate IV setups, the Morgan Lens Delivery Set may be attached to one or two Morgan Lenses.

Fluid administration tubing only—lenses not included


The Medi-Duct®

is an ocular fluid management system designed to make irrigation with the Morgan Lens even more convenient. Its super-absorbent wick carries the irrigation solution away from the patient for easy collection and disposal.

The Morgan Lens, Morgan Lens Delivery Set, Medi-Duct and MorTan's packaging materials are not made with natural rubber latex.

Symbol Glossary		
Symbol	Title Reference	Description
	Manufacturer 5.1.1[1]	Indicates the medical device manufacturer
	Authorized representative in the European Community/ European Union 5.1.2[1]	Indicates the authorized representative in the European Community/ European Union
	Date of Manufacture 5.1.3[1]	Indicates the date when the medical device was manufactured
	Use-by date 5.1.4[1]	Indicates the date after which the medical device is not to be used
	Batch code 5.1.5[1]	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalogue number 5.1.6[1]	Indicates the manufacturer's catalogue number so that the medical device can be identified
	Distributor 5.1.9[1]	Indicates the entity distributing the medical device into the locale
	Country of Manufacture 5.1.11[1]	To identify the country of manufacture of products. The date of manufacture may also appear adjacent to symbol.
	Sterilized using irradiation 5.2.4[1]	Indicates a medical device that has been sterilized using irradiation.
	Do not re-sterilize 5.2.6[1]	Indicates a medical device that is not to be re-sterilized
	Do not use if package is damaged and consult instructions for use 5.2.8[1]	Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information
	Single sterile barrier system 5.2.11[1]	Indicates a single sterile barrier system
	Single sterile barrier system with protective packaging outside 5.2.14[1]	Indicates a single sterile barrier system with protective packaging outside
	Do not re-use 5.4.2[1]	Indicates a medical device that is intended for one single use only
	Consult Instructions for Use 5.4.3[1]	Indicates the need for the user to consult the instructions for use
	Negation of - Contains or presence of natural rubber latex 5.4.5[1]; Annex B.2[1]	Indicates there IS NOT the presence of dry natural rubber or natural rubber latex as a material of construction within the medical device or the packaging of a medical device
	Medical device 5.7.7[1]	Indicates the item is a medical device
	Unique device identifier 5.7.10[1]	Indicates a carrier that contains unique device identifier information
	Prescription only 21 CFR 801.109	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	CE Marking MDR 2017/745 Article 20	Signifies European technical conformity
	Swiss authorized representative Section 3[2]	Indicates the authorized representative in Switzerland
	UK Responsible Person (No Applicable Reference)	Indicates the responsible person in the United Kingdom
	Quantity (No Applicable Reference)	Indicates the quantity of devices in the package

[1] ISO 15223-1:2021 Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements

[2] Swissmedic, Obligations Economic Operators CH, VM-ID: MU600_00_016e / V4.0 / mea / pmi / 06.04.2023

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MorTan, Inc 329 East Pine Street. Missoula, MT